



Bure™
StandTall Walker



New possibilities
for patients.
New conditions for
carers like you.

BureTM

StandTall Walker

The StandTall Walker has had a central function in healthcare even though it has lived life out of the limelight. But now GATE has thrown out all its old solutions and started again from scratch. The result is the Bure StandTall Walker with its unique ergonomics.

Improved general condition

Apart from being an excellent aid for getting patients up and about, the walker is also a valuable rehabilitation tool. Rehabilitation using the walker keeps the patient's muscles and alimentary canal active and generally contributes to improving the general condition of the otherwise bedridden. It makes it easier for infirm patients or those with weak legs to walk, and some can even manage without the support of healthcare personnel.

Improved work ergonomics

The walker creates conditions for much better working ergonomics for those who care for, train and work together with patients. For example, the open design provides increased manoeuvring space and makes access easier when getting up or dressing. The Bure StandTall Walker contributes to reduced workloads and is an investment in personnel health.

Enhanced stability

The Bure StandTall Walker has anatomically designed cushions at a slight incline that help shoulders and arms to relax. Its design helps to shift the focus from muscles in the hands, arms and shoulders to larger muscle groups such as leg and chest muscles. Stability is improved and the walker is easier to steer, at the same time as trunk bending is cushioned and posture improved.

NEW

Bure Double 2.0

EL

The new Bure Double is equipped with twin electric drives for raising and lowering the walker and widening the walking frame. The new version has a significantly lower height above the walking frame for better access and a higher user weight (150 kg) compared to the earlier version.

Sustainable

The StandTall Walker has an extra reinforcement band around the arm rest cushions to prevent them from splitting from impacts to the side. Armrests in PU foam are available as an alternative for additional comfort/hygiene.

Easily reachable hand controller

The hand controller is located such that both of the carer's and user's hands can easily reach it.

Flexible stepless adjustment

Electrically driven height adjustment makes it easy to pre-set the walker to the desired height. Once the patient is standing, the height can be finely adjusted so the walker is comfortable to use and provides the best support for walking.

Power Rise

Bure Double provides excellent assistance when helping patients stand up. The strain on personnel, and thus repetitive strain injuries, is minimized. Some patients can even get up on their own with the assistance of the walker.

Eminently suitable at home

Because Bure Double is so flexible it is eminently suitable for the home environment. With the frame arms in the inner position the StandTall Walker is narrow, easy to manoeuvre and able to negotiate most doorways.

Height-adjustable handles

See accessories.

Ergonomic handles

Ergonomically designed handles are steplessly adjustable in every direction for optimum individual comfort. This is extremely important for stroke patients who can thus reach their "usual" hand position.

Stable

The walker has anatomical cushions that provide relief for shoulders, arms and the neck, and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Open design

The open design provides the same stability but with greater room for manoeuvre. The patient does not feel shut in and it makes it easier to help with standing up and dressing.

Low height

Low ground clearance makes it easy to get close to patients e.g. when they are sitting on their beds.

Lockable castors

All four castors are individually lockable as standard on Bure Double.

Electrically driven frame widening

Electrically driven frame widening makes it easy to move the walker right up close to users without worrying about wheelchair width, chair width or bed design. It also helps personnel assist patients to stand up. With the frame in its outer position, the walker is even more stable and can be used to best advantage by the user and it also helps patients who need extra walking space.

Waterproofing

All Bure Double electrical components have IP 65 international protection ratings and can therefore be used in wet spaces.

All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors). Examples of possible parts combinations are shown below.

Technical specifications

CE marked according to MDD 93/42/EEC. | Tested and approved according to EN 11199-3/ 2005

Part No.	Product	Width	Max user weight	Length	Height	Arm rest/internal dimension
56-303/75	Bure Double 2.0 EI, 75mm castors	64-95 cm	150 kg	87 cm	85-125 cm	25-50 cm
56-303	Bure Double 2.0 EI, 125mm castors	64-95 cm	150 kg	87 cm	90-130 cm	25-50 cm

Bure Ortho/Space

ORTHO | ORTHO S

SPACE

Ortho is especially designed for orthopaedics departments and situations that demand extra legroom and room for manoeuvre for the user. Ortho has a powerful electric motor for powered height adjustment. Walker heights are thus adjustable with the patient in situ and the walker can be used as a stand-up aid. Ortho is also available in a small version – Ortho S (only available in the GAS version) – with a max frame width of only 65 cm.

In addition to Ortho's advantages, Bure Space also has extra room to get up close to wheelchairs or similar. The base frame has been made wider and the internal width extended to 67 cm. Bure Space is available with electrically operated or hydraulic height adjustment.

Ergonomic handles

Stepless adjustment in every direction.

Easily reachable hand controller

The hand controller is located such that both of the carer's and user's hands can easily reach it.

Waterproofing

All Bure Double electrical components have IP 65 international protection ratings and can therefore be used in wet spaces.

Power Rise

Bure Ortho provides excellent assistance for helping users stand up.

Unrestricted front legroom

Users with contractures and patients with full leg casts can keep their legs straight and still bring the walker close to them.

SPACE



More space

Bure Space has a base frame with internal width of 67cm - more space to get close with a wheelchair.

Height-adjustable handles

See accessories.

Low minimum height

Only 85 cm (when 75 mm castors are used).

Stable

The walker has anatomical pads that provide relief for shoulders, arms and the neck and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Flexible stepless adjustment

This makes it easy to pre-set the walker to the desired height.

Low base

This enables it to fit under most bed models.

Open design

Provides the same stability but with greater room for manoeuvres.

Brakes and directional locks

Located so that the user can reach them from the walking position

All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors). Examples of possible parts combinations are shown below.

Technical specifications

CE CE marked according to MDD 93/42/EEC. | Tested and approved according to EN 11199-3/ 2005

Part No.	Product	Width	Max user weight	Length	Height	Armrest
56-317T	Bure Ortho, 125mm casters	78 cm	150 kg	80 cm	90-130 cm	25-50 cm
56-318	Bure Ortho Gas, 125mm casters	78 cm	150 kg	80 cm	90-130 cm	25-50 cm
56-318-S	Bure Ortho S Gas, 125mm casters	65 cm	150 kg	79 cm	90-130 cm	30-45 cm
56-304	Bure Space EI, 125 mm casters	78 cm	150 kg	79 cm	90-130 cm	25-50 cm
56-307	Bure Space Gas, 125 mm casters	78 cm	150 kg	79 cm	90-130 cm	25-50 cm

Bure Standard Bure S

EL | GAS | MANUAL

Bure Standard and Bure S are available in three versions; EL | GAS | and MANUAL. Bure S is a narrower version of the Bure Standard specially adapted to work in confined spaces. Both models have stepless forearm support height adjustment.

Ergonomic handles

Ergonomically designed handles that are steplessly adjustable in every direction for optimum individual comfort. This is extremely important for stroke patients who can thus reach their "usual" hand position.

Easily reachable hand controller

The hand controller is located such that both of the carer's or user's hands can easily reach it.

Power Rise

Bure Standard EL and Bure S provide excellent assistance when helping patients and users stand. The strain on personnel, and thus repetitive strain injuries, are minimized. Some patients can even get up on their own with the assistance of the StandTall Walker.

Eminently suitable at home

Bure S is perfectly suited to confined spaces such as at home or in some care accommodation. Bure S will pass through most doorways that are 70 cm or wider.

BURE S



EL

Bure Standard and Bure S EL are ideal for users who need a walkers equipped with power rise and which are capable of height adjustment once a patient is standing.

GAS

Bure Standard and Bure S GAS are ideal for users who need a quality walker capable of stepless height adjustment but who do not need assistance from a power rise actuator.



MANUELL

Bure Standard and Bure S MANUELL are great for users who need good, individually adaptable walkers but who do not need to adjust the height very often and who have no need of assistance from a power rise actuator.



Height-adjustable handles
See accessories.

Sustainable
StandTall Walkers are designed to last. For example, they have extra reinforcement bands around the arm rest cushions to prevent them from splitting from impacts to the side.

Stable
The walker has anatomical pads that provide relief for shoulders, arms and the neck and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Flexible stepless adjustment
Electrically driven height adjustment makes it easy to pre-set the StandTall Walker to the desired height. Once the patient is standing in position the height can be steplessly adjusted for comfort and the proper support for an optimal walking configuration.

Brakes and directional locks
There are optional one or two-handed brakes that can be positioned so that the user can reach them from the walking position. The handbrake acts on the rear castors and the brake assembly also includes directional locks that make it easier for the user to stay on course when walking.

Open design
The open design provides the same stability but with greater room for manoeuvre. The patient does not feel shut in and it makes assistance with standing up and dressing easier.

Waterproofing
All Bure Standard EL and Bure S EL parts have IP 65 international protection ratings and can therefore be used in wet spaces.

All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors). Examples of possible parts combinations are shown below.

Technical specifications

CE marked according to MDD 93/42/EEC. | Tested and approved according to EN 11199-3/ 2005

Part No.	Product	Width	Max user weight	Length	Height	Arm rest/internal dimension
56-312T	Bure Standard EL	78 cm	150 kg	80 cm	90-130 cm	25-50 cm
56-311	Bure Standard GAS	78 cm	150 kg	80 cm	90-129 cm	25-50 cm
56-310	Bure Standard MANUELL	78 cm	150 kg	80 cm	90-133 cm	25-50 cm
56-302T	Bure S EL	66 cm	150 kg	79 cm	90-130 cm	30-45 cm
56-301	Bure S GAS	66 cm	150 kg	79 cm	90-129 cm	30-45 cm
56-300	Bure S MANUAL	66 cm	150 kg	79 cm	90-133 cm	30-45 cm

Bure Low

EL | GAS

Bure is now also available in a lower model for short users. The model is the same as Bure Std/S except that the min/max heights are significantly lower (70-95cm with 75mm casters). Bure Low is available with two different lower frames, Std and S. Both versions are available with electric height adjustment or height adjustment using gas struts.

Ergonomic handles

Ergonomically designed handles that are steplessly adjustable in every direction for optimum individual adaptation. This is extremely important for stroke patients who can thus reach their "usual" hand position.

Easily reachable hand controls

The hand controls are located such that both of the carer's or user's hands can easily reach them.

Power Rise

Bure Standard EL and Bure S provide excellent assistance when helping patients and users stand. Personnel burdens, and thus repetitive strain injuries, are minimized. Some patients can even get up on their own with the assistance of the Standtall Walker.

Eminently suitable at home

Bure S is perfectly suited to confined spaces such as at home or in some care accommodation. Bure S will pass through most doorways that are 70 cm or wider.

Height-adjustable handles

See accessories.

Sustainable

The walker is designed to last a long while. For example it has extra reinforcement bands around the arm rest cushions to prevent them from splitting from impacts to the side.

Stable

The walker has anatomical pads that provide relief for shoulders, arms and the neck and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Flexible stepless adjustment

Electrically driven height adjustment makes it easy to pre-set the Standtall Walker to the desired height. Once the patient is standing in position the height can be steplessly adjusted for comfort and the proper support for an optimal walking configuration.

Brakes and directional locks


There are one or two optional hand brakes that can be positioned so that the user can reach them from the walking position. The handbrake acts on the rear castors and the brake assembly also includes directional locks that make it easier for the user to stay on course when walking.

Open design

The open design provides the same stability but with greater room for manoeuvre. The patient does not feel shut in and it makes assistance with standing up and dressing easier.

Waterproofing

All Bure Standard EL and Bure S EL parts have IP 65 international protection ratings and can therefore be used in wet spaces.



Bure Std

Bure Low



All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors). Examples of possible parts combinations are shown below.

Technical specifications

CE CE marked according to MDD 93/42/EEC. | Tested and approved according to EN 11199-3/ 2005

Part No.	Product	Width	Max user weight	Length	Height	Arm rest (internal dimension)
56-321/75	Bure Low EL, 75 mm casters	78 cm	150 kg	80 cm	70-95 cm	25-50 cm
56-321/75-S	Bure Low S, 75 mm casters	66 cm	150 kg	79 cm	70-95 cm	30-45 cm
56-322/75	Bure Low Gas, 75 mm casters	78 cm	150 kg	80 cm	70-95 cm	25-50 cm
56-322/75-S	Bure Low S Gas, 75 mm casters	66 cm	150 kg	79 cm	70-95 cm	30-45 cm

Bure XL

EL

Bure XL was designed for the bigger, heavier patient and therefore provides extra space and strength. Because extra powerful aids are required in the care and nursing of heavier patients, Bure XL is always supplied with electric drives.

Ergonomic handles

Ergonomically designed handles that are steplessly adjustable in every direction for optimum individual comfort. This is extremely important for stroke patients who can thus reach their "usual" hand position.

Easily reachable hand controller

The hand controller is located such that both of the carer's or user's hands can easily reach it.

Power Rise

Bure XL is a good aid for helping patients and users stand and reduces strain on employees. Some patients can even get up on their own with the aid of the walker.

Height-adjustable handles

See accessories.

Sustainable

StandTall Walkers are designed to last. For example, they have extra reinforcement bands around the arm rest cushions to prevent them from splitting from impacts to the side.

Stable

The walker has anatomical cushions that provide relief for shoulders, arms and the neck, and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Flexible stepless adjustment

Electrically driven height adjustment makes it easy to pre-set the walker to the desired height. Once the patient is standing in position the height can be steplessly adjusted for comfort and the proper support for an optimal walking configuration.

Open design

Provides the same stability but with greater room for manoeuvres.

Waterproofing

All Bure XL electrical components have IP 65 international protection ratings and can therefore be used in wet spaces.

Hospitals and nursing homes

Because of its size, Bure XL is best suited to hospitals and nursing homes.

Designed for weight

The Bure XL is designed to provide support for heavier patients up to 240 kg. This StandTall Walker has also been widened to provide good stability and plenty of space. The inner dimension is 44–62 cm. The extra leg room makes it especially easy for patients who need extra walking space due to e.g. leg spasms.

All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors). Examples of possible parts combinations are shown below.

Technical specifications

CE CE marked according to MDD 93/42/EEC. | Tested and approved according to EN 11199-3/ 2005

Part No.	Product	Width	Max user weight	Length	Height	Arm rest/internal dimension
56-315	Bure XL EL	88 cm	240 kg	78 cm	91-131 cm	44-62 cm

Accessories

NEW FEATURES

New, improved platform

We have developed a new platform to replace the current aluminium model (does not apply to the Bure Ortho or the Bure XL platform). The new platform is hung from the front of the walker and has its own rear castors to facilitate manoeuvring. The new platform is not only lighter and smaller, it is also simple and easy to use.



Handlebar for easier manoeuvring

We've developed a new handlebar to make it easier for carers to help wheelchair users in the best, most convenient way. The handlebar makes it easy to manoeuvre our walker, especially when the user is standing on the platform.

The handlebar is easy to fit and can be retrofitted on all Bure Walkers (except for Bure XL).



Shin supports

A shin support is now available to further extend and improve the use of our standard walkers for help with standing and mobility. The accessory can be used with all walkers except Rise&Go (DB) and Bure XL.



Because of the tip risk depending on the patient and walker model, we recommend the accessory only be used in combination with a platform and/or stirrups/tip protection.

GATE PLUS

Did you know GATE now offers additional functions on every electric walker? The new functions are integrated into the hand controller/PLUS (accessory). We call these new functions GATE PLUS.

Hand controller/PLUS is suitable for all walkers (not Bure XL) produced after 25 November 2015.

Hand controller/PLUS can be specified when placing a new order or ordered as an accessory and will provide the following functions:

- Multiple speeds – there is a fast mode in addition to the current standard speed (+40% compared to standard speed).
- Battery indicator in the hand controller – makes it easy to see when it's time for a recharge.
- Service indicator – helps you keep track of when your walker needs service and maintenance (at least once a year).



Hygiene Management – new armrests pads in PU

As an alternative to our current cold foam armrest pads, we are launching pads in moulded polyurethane. Their primary advantages are the lack of perforations or edges that allow water and bacteria to penetrate

and remain in the surface after cleaning and/or disinfecting. This makes it easier to keep the walker clean between each user.

Our new polyurethane armrest pads are available for all walkers. Polyurethane shin supports are also available for Rise & Go (DB). The new pads can be specified when placing new orders or ordered as retrofit accessories. Pls note that the armrest pads PU will add 5cm on the overall width compared to measures given.



Walker movement controlled by buttons in the handles.

We now offer remote control also through buttons in the handles.

We offer push button or rocker switch control. All functions can be operated in this way. The accessory is also retrofittable and when necessary it can be used in combination with the regular hand controller.



Gait training belt STD

Many have asked for a functional Gait training belt that can be used with our standard range of StandTall walkers. Our new Gait training belt is attached to 4 points providing good support for training and safety in case of falling.

Max userweight = 150 kg.

Accessories

Drip holder with attachment

Provides great freedom for the patient to move around at home or on the ward. Supplied with attachment fitting.

Hand controller

GATE PLUS provides a number of extra functions such as multiple speeds and or service/battery indicator in the hand controller.

PU foam armrest pads

PU foam armrest pads can be specified with the order or selected retroactively. The pads have no perforations or edges that accumulate dirt. An excellent alternative where high standards of hygiene management are demanded.

Pls note that the armrest pads PU will add 5cm on the overall width compared to measures given.

Soft basket with attachment

Basket with attachment has room for all the important stuff.



Stirrup

Acts as a counter weight and provides extra stability and safety when the walker is used for stand up assistance. Standard på Rise & GO. Not for use on walkers with lower frame widening.

Hand brake kit, one or two brakes

Adds a hand brake to the StandTall Walker so that the user is able to halt the walker's movement during walking exercises. Castors with directional locks included for extra safety/assistance. The one-handed brake is available left and right-handed.

When the handbrake kit is installed, base frame width is increased by +2cm.



Control for Raise/Lower and handle widening

We offer push button or rocker switch control. All functions can be operated in this way. The accessory is also retrofittable and when necessary it can be used in combination with the regular hand controller.

Adjustable side supports

Prevents arms from sliding from arm rest cushions. Steplessly adjustable along the StandTall Walker's sides to further stabilize walking. The side support is also available in a PU version.

Oxygen bottle holder with bracket

Means patients are no longer confined to bed.

Shin supports

Leg supports make standing and mobility easier. However, they should be used with a platform and stirrup or tip protection due to the risk of tipping backwards.

Platform

The platform is an excellent addition when you want to use the walker for standing and mobility. Available in a short or long version.

Heel strap


Hand brake kit

Directionally lockable castor

125 mm
100 mm
75 mm

Castor with brake

125 mm
100 mm
75 mm



Accessories



Drip rod complete with attachment



Side support/Side support PU



Platform short/long



Platform, short/long, castors (SW)



Handbrake kit DB



Handbrake kit standard



Castor with brake 75, 100, 125 mm



Castor with directional locking 75, 100, 125 mm



Shin supports



Soft basket



Tube holder incl bracket



Armrest cushions, PU



Handlebar



Heel strap to platform



Control/handle
- rocker switch




Control/handle
- push button switch



Tip protection - backwards



Stirrup



Gait training belt STD

Did you know that...

HINT



The new Bure DB now has lower ground clearance, which makes it significantly easier to fit under beds and thus get closer to the user. When fitted with 75 mm castors, Bure DB is especially low.



If you choose a walker with electrically adjustable widening you are also choosing increased flexibility. You can get close to the patient regardless of wheelchair or chair width. Also note that in most cases the platform will fit between a wheelchair's front castors.



The platform is an excellent addition when you want to use the walker for mobility, standing up and patient transits. The platform is narrow enough for use on most wheelchairs.

Quick list

ACCESSORIES

Designation	Part #	Bure Std/S	Bure Low	Bure Ortho	Bure Space	Bure DB 2.0	Bure XL
Armrest cushions PU, Std/Low/DB/Ortho/Space/XL	57-361-PAR	X	X	X	X	X	X
Armrest cushions (Std) Foam, Std/Low/DB/Ortho/Space	56-366-PAR	X	X	X	X	X	X
Drop rod complete with holder	56-338	X	X	X	X	X	X
Electric motor kit + battery incl battery shelf, Low	56-332T-LOW		X				
Electric motor kit + battery incl battery shelf, Std	56-332T	X					
Gas strut kit incl complete, Low	56-330-LOW		X				
Gas strut kit incl complete, Std	56-330	X					
Gait training belt STD	56-388	X	X	X	X	X	
Castor							
Castor with brake 75mm	56-387	X	X	X	X	X	X
Castor with brake 100mm	56-358	X	X	X	X	X	X
Castor with brake 125mm	56-361	X	X	X	X	X	X
Castor with directional locking 75mm	57-043	X	X	X	X	X	X
Castor with directional locking 100mm	56-355	X	X	X	X	X	X
Castor with directional locking 125mm	56-359	X	X	X	X	X	X
Heel strap/Platform	58-328-1	X	X	X	X	X	
Shin support	56-352	X	X	X	X	X	
Handlebar	57-368	X	X	X	X	X	
Manual brake							
Manual brake, 2 hands, Std/Low/Ortho/Space/XL	56-339-1-0	X	X	X	X		X
Manual brake, 1 hands, Std/Low/Ortho/Space/XL	56-339-1-1	X	X	X	X		X
Manual brake, 2 hands, DB 2.0	56-339-2-0					X	
Manual brake, 1 hands, DB 2.0	56-339-2-1					X	
Soft basket incl bracket	56-389	X	X			X	X
PLUS							
Hand controller PLUS incl holder, Std/Low/Ortho/Space	58-375-1	X	X	X	X		
Hand controller PLUS incl holder, DB/Rise & Go DB	58-375-2					X	
Side support	56-337	X	X	X	X	X	X
Side support, PU	56-337-PU	X	X	X	X	X	X
Steering/Handle							
Button kit FK, Up/Down, Normal speed, one htg*	58-387-1	X	X	X	X	X	
Button kit FK, Up/Down, High speed, one htg*	58-387-2	X	X	X	X	X	
Button kit FK, Widening, one htg*	58-387-3					X	
Button kit TK, Up/Down, Normal speed, two htg*	58-388-1	X	X	X	X	X	
Button kit TK, Up/Down, High speed, two htg*	58-388-2	X	X	X	X	X	
Button kit TK, Widening, two htg*	58-388-3					X	
Y cable for button kit – one per walker*	57-022-1	X	X	X	X	X	
Y cable for use with hand controller together with button kit	57-022-3	X	X	X	X	X	
Platforms***							
Platform short, STD	56-333-STD, 56-333/SW-STD	X	X				
Platform short, DB	56-333-DB, 56-333/SW-DB					X	
Platform short, SPACE	56-333-SPACE, 56-333/SW-SPACE				X		
Platform long, STD	56-383-STD, 56-383/SW-STD	X	X				
Platform long, DB	56-383-DB, 56-383/SW-DB					X	
Platform long, SPACE	56-383-SPACE, 56-333/SW-SPACE				X		
Platform XL	56-386						X
Platform Aluminium short, Ortho	56-392			X			
Platform Aluminium long, Ortho	56-394			X			
Tip protection – backwards, Std/Low/Ortho/Space**	57-387	X	X	X	X		
Stirrup, Std/Rise & Go	56-331-1	X	X				
Tube holder incl bracket	56-336	X	X	X	X	X	X

*Tools required for drilling holes in frame when retrofitting.


**Can not be used with Manual brake.

*** /SW is variant with castors.

BureTM

StandTall Walker

- User- friendly design with smooth lines that do not confine the patient.
- Structure improves stability and control.
- Available in eight basic models with great possibilities for individual customization.
- A range of functional accessories complement individual fit perfectly.
- Design and functions that create great ergonomics and make day-to-day work easier.
- Electrically driven power rise assistance means that patients can often actively participate in standing up.
- Electrically driven frame widening creates greater freedom and improves the standing up and patient transit functions.



WE DO BETTER.

Want to know more?

Contact us: tel +46 (0)371-318 00
Gate Rehab Development AB
Industrigatan 2 | SE-333 75 Reftale SWEDEN
info@gaterd.com | www.gaterd.com



GATE COLLECTION