

# Amfibi Double

EL

Amfibi Double is GATE's most advanced and versatile hygiene chair. It is equipped with twin electric drives, electrically operated raising and lowering and electrically operated tilt. The chair can be adjusted steplessly even when the patient is in the chair. Together with the soft, relaxing relief offered by the seat and back, the chair's functions provide superlative comfort for the patient and a minimum of strain for personnel.

**Best possible comfort**

## Tilt function

Amfibi Double can be inclined 5 degrees forward to 25 degrees rearward. This makes it easy to use for personal hygiene and when sitting down or getting up. The tilt function also improves the working environment and reduces strain on personnel.

## Ergonomic seat

The Amfibi Double's seat is adapted for an upright posture that does not close the pelvis. It minimizes pressure and is designed to suit both men and women. The seat is easy to remove and clean, while its rounded front edge contributes to a comfortable leg position.

## Plastic or metal

A guide is available for Amfibi as an accessory. Note that the guide is available in different versions to suit plastic or metal bedpans. See accessories.

## Swing-away, detachable footrests

Flip-up, swing-away, detachable footrests with heel restraints to prevent the patient's feet from sliding off in the tilt position are fitted as standard.

## Adjustable seat height

An infinitely variable seat height from 53 all the way to 82 cm means the chair will fit every patient and reduce staff workload to a minimum.

## Individual fit

Amfibi Double has a neck support that is adjustable in height and depth for comfortable support and individual fit.

## Easy control

Simply press the remote up/down and tilt buttons to operate. The cabled remote is sealed to a high protection class (IP65).

## Adjustable back-strap

The soft, comfortable back-strap relieves pressure on bony areas. It is easily adjustable for best support and seat depth. Machine washable.

## Lockable armrests

The armrests can be flipped up and locked in position to make it easier for the patient to assist with the transfer.

## Ergonomics

The chair liberates the user from effort completely, while no manual exertion is required on the part of the carer during any of the manoeuvres.

## Durable materials

Manufactured from corrosion-proof materials.

## Lockable castors

All four castors are lockable.



For accessories see page 7.

## Technical specifications

CE marked according to MDD 93/42/EEC.

Part No.	Product	Width, external	Max user weight	Seat height	Width between armrests	Seat depth
56-203DB	Amfibi Double	62 cm	150 kg	53-82 cm	52 cm	45 cm